



AURA & SPECTRA Light Engines®

Microscope adapter options

| Manufacturer | Model | Lumencor Light Engine Adapter | |
|-------------------------------|------------------------------|-------------------------------|---|
| | | AURA part number | SPECTRA part number |
| Carl Zeiss | AxioExaminer | 82-10001 | 82-10009 |
| | Axiolmager | | |
| | AxioObserver | | |
| | AxioScope A1 | | |
| | Axiophot | | |
| | Axioplan 1, 2 | | |
| | Axioskop 2, 2FS, 40 | | |
| | AxioTech | | |
| | Axiovert 200 | | |
| Leica | DM2000, DM6000, DMI4000 | Please inquire | 82-10016 |
| | DMIL | | |
| | DMIRB, DMIRM, DMIRE2 | | |
| | DMLB | | |
| | DMLFS | | |
| | DMR, DMRXA, DMRXA2 | | |
| Nikon | 50i, 80i, 90i | 82-10003 | 82-10011 |
| | E400, E600, E800 | | |
| | FN-1 | | |
| | SMZ1500, SMZ800 | | |
| | TE2000, all other "TE" | | |
| | Ti | | |
| | TS100 | | |
| | AX70 | | |
| Olympus | BX series: BX50, BX60 | 82-10002 | 82-10010 |
| | BX2 series: BX41, BX51, BX61 | | |
| | BX3 series: BX43, BX46, BX53 | | |
| | CKX41 | | |
| | CK40 | | |
| Liquid light guide | 3 mm core Lumatec | 82-10006 | 82-10012 |
| | 5 mm core Lumatec | 82-10020 | 82-10022 |
| | custom | upon request | upon request |
| 1 mm quartz fiber | SMA or FC/PC connector | 82-10004 | 82-10015 |
| Yokogawa Confocal Head | CSU 10 Spinning Disk | 82-10017 | Reverse dichroics for GFP&mCherry and CFP&YFP fluorophor sets available Please inquire |

Custom adapters available upon request, (e.g. fiber bundle, projection optics), please inquire